

EEG → Brain wave medical test

- ↳ Electrodes via scalp
- ↳ Hair is dry & oil free. Since many electrodes inserted

Normal Brain wave → 30 Hz

Epilepsy Brain → Burst of abnormal spikes

EEG of babies is irregular → outlier maybe

30-60 min to perform EEG

Not good for a person with epilepsy.

Tumor or stroke → very low EEG freq

helps give a general overview but does not specify or pinpoint the exact location.

16 - 25 electrodes. while testing





- Find area to put the electrodes
- Forehead-Nose maybe.
- Electrode put on chest to record the electrogram.

### Sites where electrodes are inserted

- F: Frontal
- FP: Frontopolar
- P: Parietal
- A: Auricular (ear electrode)
- T: Temporal
- C: Central
- O: Occipital

→ we can even more localise the location

Such as T3, T4, T5

- Even number → Right side of head.
- Odd Number → Left side of head.

### Ways to connect EEG

- Records brain waves using amplifiers.
- on a bipolar montage → electrode input 1 of one channel & input 2 of another

→ Being able to tell the exact location of electrical activity ("Localisation") is critical, being able to interpret EEG meaningfully.

### → Localising in Bipolar montages

- By identifying the phase reversal.
- A deflection within a chain pointing to off direction.

### → In referential montage


- All waves move in a positive/same direction
- when +ve wave and we see a change then deflection downwards
- Electrode with maximum negative activity in referential montage



## HOW TEST IS SET UP

- ↳ Relative position lines to pin point exact locations.
- ↳ 21 recording electrodes + 1 ground electrodes
- And one EKG electrode for eye movement

**Amplification** → Measure activity of one electrode relative to another.



Common data is eliminated and hence makes easier.

→ Improves signal to noise ratio in recording

Filtering is done to just get 1-30 Hz range.

**Filtering**

- Low freq (Slow waves)
- High freq (Fast)
- Notch filter (Narrow region)

Identifying the lead region with phase reverse.

→ Problem → Phase cancellation of biological activity → might just cancel each other  
• Leads to false localisation


Detects a large group of neurons that are active at the same time.

→ Measures Post Synaptic Potentials (Not action potential)




## DIPOLE PRINCIPLE

diff charged electrodes at a distance




↓  
This graph is due to  
diff in charge b/w these  
two areas

## Excitatory Post Synaptic Potential




Inhibitory Post Synaptic Potential  
↳ what if the dipoles we created is  
inverted



These 2  
create  
differences  
in charge

→ we can measure this difference but  
it is too small of a number to work  
on

→ clusters of pyramidal neurons are  
stimulated and create an overall  
change -






EEG Has 32 electrodes like this in multiple regions.

→ Need for pairing up of electrodes

→ Need to find perfect pair combo.


① perspective of one electrode to all other electrodes



compare every electrode to same electrode

01-CZ		} common self montay
P4-CZ		
P3-CZ		
F4-CZ		
F5-CZ		
T8-CZ		

② Other way is to compare every electrode relative to its adjacent



### 5 Key Bands of freq signals

Sleep → Delta →  $< 4$

Relaxing → Theta →  $4-8$

Attention → Alpha →  $8-13$

Active Mind → Beta →  $13-35$

High Concent → Gamma →  $> 35$

Interictal Epileptiform Discharges  
↳ consists of sharp spike waves

Normal  $\rightarrow$  20 - 70 ms

Spike  $\rightarrow$  70 - 200 ms.

Artifacts  $\rightarrow$  Non brain signal that  
contribute to the EEG.

- $\rightarrow$  Eye Blink
- $\rightarrow$  muscle movements
- $\rightarrow$  Power line interference